



ADD: NO. 5 XINGFU RD, ZONE A, SHISHAN INDUSTRIAL PARK, NANHAI DISTRICT, FOSHAN CITY, GUANGDONG, P.R. CHINA.

E-mail: fsdongfang@fsdfmed.com

Http: //www.fsdfmed.com

PC: 528225

TEL: +86-757-81202183 81202112 88368396 81202181

FAX: +86-757-81202189



SINCE 1990

Introduction

ESTABLISHED IN 1990, FOSHAN DONGFANG MEDICAL EQUIPMENT MANUFACTORY (LTD) IS AN ISO9001-2015 AND ISO13485:2012 CERTIFIED FACTORY WHO MANUFACTURES VARIOUS KINDS OF REHABILITATION EQUIPMENTS, INCLUDING WHEELCHAIRS, WALKING AIDS AIDS AND HOSPITAL BEDS, ETC.

COVERING AN AREA OF 30,000 SQUARE METERS, THE FACTORY HAS OVER 700 EMPLOYEES, PRODUCES OVER 500,000 WHEELCHAIRS, 2,000,000 WALKERS, AND 18,000 HOSPITAL BEDS PER YEAR; AND IS EXPORTING TO OVER 80 COUNTRIES WORLDWIDE.

THE FACTORY CONCENTRATES ITSELF ON THE PRODUCT QUALITY AND INNOVATION. IN 1999, IT WAS THE FIRST DOMESTIC WHEELCHAIR FACTORY THAT OBTAINED THE CE QUALITY APPROVAL. IN 2003, IT GOT THE 510K REGISTRATION FROM USA FDA. IN 2007, IT PASSED "SG" QUALITY APPROVAL FROM CPSA (CONSUMERPRODUCTS SAFETY ASSOCIATION OF JAPAN).

AFTER 30 YEARS OF HARD WORKING, FOSHAN DONGFANG MEDICAL HAS EARNED A HIGH REPUTATION FOR ITS EXCELLENT QUALITY AND SERVICE IN GLOBAL MARKET. AND IT ALWAYS ENDEAVORS TO FURTHER IMPROVE PRODUCT QUALITY IN ORDER TO BETTER USERS' LIFE STANDARDS.



Certificates



Function Code

FS253LCHQ-46

(-46): SEAT WIDTH 46cm(18")
Q: QUICK RELEASE AXLE
H: DRUM BRAKE
C: ELEVATING LEGREST
L: ALUMINUM
MODEL NO.



A: DROP BACK HANDLE B: MAG WHEEL C: ELEVATING LEGREST E: DETACHABLE PUSH BAR



G: RECLINING HIGH BACK H: DRUM BRAKE J: UNITED BRAKE L: ALUMINUM P: PU TYRE



Q: QUICK RELEASE AXLE S: STAINLESS STEEL U: U-SHAPE COMMODE SEAT X: INCLINED ARMREST Y: HARD CUSHION SEAT

FS3220WMF2

F2: DETACHABLE GUARDRAILS
M: WOODEN HEAD AND FOOT BED-BOARD
W: PEDAL LOCKING CASTORS
MODEL NO.



B: PLASTIC BED BOARD F1: DETACHABLE LONG GUARDRAILS F2: DETACHABLE GUARDRAILS F3: FOLDABLE LONG GUARDRAILS F4: SECTIONAL TYPE ELEVATING GUARDRAILS



G: FLAT METAL BED BOARD M: WOODEN HEAD AND FOOT BED-BOARD W: PEDAL LOCKING CASTORS Z: CENTRAL LOCKING CASTORS

Contents

59-74

Commode Series

01-20

Aluminum Wheelchair Series

75-80

Bath Bench Series

21-26

Reclining Wheelchair Series

81-102

Walking Aids Series

27-38

Power Wheelchair Series

103-118

Hospital Bed Series

39-46

Leisure & Sports Wheelchair Series

47-58

Steel Wheelchair Series

119-128

Rehabilitation Aid Series

129-132

Spare Parts/Colour Sample

Aluminum Wheelchair Series



FS205LHQ

- › Aluminum chair frame
- › Detachable footrest
- › Solid castor
- › Locking drum brake
- › Anti-tippers
- › Carrying wheel
- › Quick release pneumatic rear wheel
- › Angle-adjustable footplate
- › Flip-up and height-adjustable inclined armrest

MAIN PARAMETERS (cm)

105	91	68	37.5	46	100
61	20	54	77x33x93	N.W.17.8kg	G.W.20.1kg

Always care for your health



FS250LCPQ

- › Aluminum chair frame
- › Flip-up and height-adjustable inclined armrest
- › Detachable elevating footrest
- › PU castors
- › Quick release pu rear wheel
- › Anti-tippers.

MAIN PARAMETERS (cm)

110	60	67	38	46	100
60	20	52	83x31x93	N.W.17.4kg	G.W.20.6kg



FS251LHPQ

- › Aluminum chair frame
- › Detachable footrest
- › PU castor
- › Quick release pu rear wheel
- › Anti-tippers
- › Locking drum brake
- › Angle-adjustable footplate
- › Flip-up and height-adjustable inclined armrest

MAIN PARAMETERS (cm)

101	90	68	39	46	100
61	20	53	80x33x92	N.W.18.1kg	G.W.20.4kg

Foshan Wheelchair



FS253LACHQ

- › Aluminum frame
- › Drop back handle
- › Locking drum brake
- › Flip-up and height-adjustable inclined armrest
- › Detachable elevating footrest
- › Quick release pneumatic rear wheel
- › Anti-tippers
- › Angle-adjustable footplate

MAIN PARAMETERS (cm)

114	94	68	42	46	100
60	20	52	85x35x78	N.W.19.3kg	G.W.21.5kg

Aluminum Wheelchair Series



FS950LBPQ

- › Aluminum chair frame
- › Detachable footrest
- › Solid castor
- › Anti-tippers
- › Quick release pu mag rear wheel
- › Rear wheel can be adjusted horizontally and vertically
- › Angle-adjustable footplate
- › Flip-up and height adjustable armrest

MAIN PARAMETERS (cm)

宽度 (Width)	深度 (Depth)	高度 (Height)	扶手高度 (Armrest Height)	脚踏板高度 (Footrest Height)	总高度 (Total Height)
98	90	66	38	46	100
61	20	52	80x34x93	N.W.16.8kg G.W.19kg	

Always care for your health



FS955L

- › Aluminum chair frame
- › Flip-up inclined armrest
- › Detachable footrest
- › Solid castor
- › Pneumatic rear wheel
- › Angle-adjustable footplate

MAIN PARAMETERS (cm)

宽度 (Width)	深度 (Depth)	高度 (Height)	扶手高度 (Armrest Height)	脚踏板高度 (Footrest Height)	总高度 (Total Height)
107	91.5	63	38	46	100
61	20	53	83x28x94	N.W.15.3kg G.W.17.7kg	



FS959LQ

- › Aluminum chair frame
- › Detachable footrest
- › PU castor
- › Quick release pneumatic rear wheel
- › Anti-tippers and carrying wheel
- › Angle-adjustable footplate
- › Flip-up and height-adjustable inclined armrest

MAIN PARAMETERS (cm)

宽度 (Width)	深度 (Depth)	高度 (Height)	扶手高度 (Armrest Height)	脚踏板高度 (Footrest Height)	总高度 (Total Height)
105	90	64	39	46	100
61	20	52.5	83x29x90	N.W.15.7kg G.W.18.3kg	

Foshan Wheelchair



FS903LQ

- › Aluminum chair frame
- › Flip-up desk armrest
- › Detachable footrest
- › Solid castor
- › Quick release pneumatic rear wheel
- › Angle-adjustable footplate

MAIN PARAMETERS (cm)

宽度 (Width)	深度 (Depth)	高度 (Height)	扶手高度 (Armrest Height)	脚踏板高度 (Footrest Height)	总高度 (Total Height)
100	90	66	38	46	100
61	20	52	80x28x91	N.W.13.9kg G.W.15.8kg	

Aluminum Wheelchair Series



AU018LPQ

- › Aluminum chair frame
- › Flip-up detachable armrest
- › Height adjustable armrest pad
- › Detachable footrest
- › PU castor,
- › Quick release PU rear wheel
- › With anti tippers

MAIN PARAMETERS (cm)

105	83	65	44	46	100
60	20	51	83x32x91	N.W.17kg	G.W.19kg

Always care for your health



FS957LQ

- › Aluminum chair frame
- › Detachable footrest
- › Solid castor
- › Quick release pneumatic rear wheel
- › With anti-tippers
- › Angle-adjustable footplate
- › Flip-up and height-adjustable inclined armrest

MAIN PARAMETERS (cm)

104	91.5	65	39	46	100
61	20	53	83x28x94	N.W.16.9kg	G.W.19.1kg



AU255LPQ

- › Aluminum chair frame
- › Flip-up and detachable armrest
- › Height adjustable armrest pad
- › Detachable footrest
- › Adjustable angle footplate
- › PU castor
- › Quick release PU rear wheel
- › With anti tipper

MAIN PARAMETERS (cm)

103	95	66	38	46	100
60	15	57	86x35x97	N.W.19.6kg	G.W.23kg

Foshan Wheelchair



FS952LCQ

- › Aluminum frame
- › Flip-up and height adjustable armrest
- › Quick release pneumatic rear wheel
- › Elevating footrest
- › Solid castor
- › Anti-tippers

MAIN PARAMETERS (cm)

109	91	67	38	46	100
61	20	52	80x34x93	N.W.17.3kg	G.W.20.1kg

Aluminum Wheelchair Series

• Foshan Wheelchair



AU974LAH-46

Aluminum Chair Frame, Fixed Armrest, Height Adjustable Armrest Pad, Detachable Footrest, Drop Back Handle, Solid Castor, Pneumatic Rear Wheel, With Drum Brake And Anti-tippers.



FS903LAJPQF9

Aluminum Chair Frame, Drop-back Handle, Flip-up Armrest, With Carrying Wheel, Detachable Footrest, With United Brake, Solid Castor, Quick Release PU Rear Wheel,



FS908LJ

Aluminum Chair Frame, Flip-up Desk Armrest, Detachable Footrest, Solid Castor, Pneumatic Rear Wheel, With United Brake



FS908LJPF25

Aluminum Chair Frame, Flip-up Footrest, Flip-up Armrest, PU Rear Wheel, Solid Castor, United Brake



FS908LAJPF21

Aluminum Chair Frame, Drop Back Handle, Flip Up Armrest, Detachable Footrest, PU Rear Wheel, With United Brake,



FS908LAJPF24

Aluminum Chair Frame, Height-adjustable Footplate, Flip-up And Height-adjustable Inclined Armrest, PU Rear Wheel, Solid Castor, Drop Back Handle, United Brake



FS208LABP

Aluminum Chair Frame, Flip-up Armrest, Detachable Footrest, Drop Back Handle, PU Mag Rear Wheel,



AU874LAHP

Aluminum Chair Frame, Fixed Armrest, Fixed Footrest, Solid Castor, PU Rear Wheel, Drop Back Handle With Drum Brake,

MAIN PARAMETERS (cm)	MODEL NO.	Dimensions and Weight													
		Width (W)	Depth (D)	Height (H)	Seat Width (SW)	Seat Depth (SD)	Backrest Width (BW)	Backrest Height (BH)	Armrest Height (AH)	Footrest Height (FH)	Footrest Width (FW)	Wheelchair Weight (kg)	DIM	NW	GW
	AU974LAH-46	104	64	46	58	15	53	89	37	100	80x34x70	17.8	19.3		
	FS903LAJPQF9	104	67	46	56	15	47	92	44	100	80x36x73	14.7	16.7		
	FS908LAJPF21	101	63	46	60	20	48.5	88.5	41	100	80.5x30x78	13.8	16.1		
	FS908LAJPF24	95	68	46	56	17.5	49	88	39	100	90x30x73	14	16		

MAIN PARAMETERS (cm)	MODEL NO.	Dimensions and Weight													
		Width (W)	Depth (D)	Height (H)	Seat Width (SW)	Seat Depth (SD)	Backrest Width (BW)	Backrest Height (BH)	Armrest Height (AH)	Footrest Height (FH)	Footrest Width (FW)	Wheelchair Weight (kg)	DIM	NW	GW
	FS908LJ	100	63	46	61	20	51	89	38	100	80x28x91	13.3	16		
	FS908LJPF25	103	64	46	56	20	52.5	91.5	41	100	81x30x96	13.3			
	FS208LABP	104	68	46	56	15	51	91	41.5	100	80x28x78	15	17.2		
	AU874LAHP	104	71	46	60	17	52	91.5	37	100	94x34x78	14.6	16		

Aluminum Wheelchair Series



FS864LJPF1

- > Aluminum chair frame
- > Fixed armrest and footrest
- > Solid castor
- > PU rear wheel
- > With united brake

MAIN PARAMETERS (cm)

102	89	66	40	46	100
61	20	48	94×28×92	N.W.12.7kg	G.W.15kg

Always care for your health



FS869LAJPF4

- > Aluminum light weight chair frame
- > Fixed inclined armrest
- > Fixed footrest
- > Solid castor
- > PU rear wheel
- > With united brake
- > Drop back handle

MAIN PARAMETERS (cm)

102	92	64	41.5	46	100
61	20	51	101×28×80	N.W.13kg	G.W.15.2kg

09

10



FS874LJPF7

- > Aluminum chair frame
- > Fixed armrest
- > Fixed footrest
- > Solid castor
- > PU rear wheel
- > United brake

MAIN PARAMETERS (cm)

94	88	61	37	46	100
56	17.5	49	90×27×91	N.W.11.8kg	G.W.14.1kg

Foshan Wheelchair



FS874LJPF5

- > Aluminum light weight chair frame
- > Fixed armrest
- > Fixed footrest
- > Solid castor
- > PU rear wheel
- > With united brake

MAIN PARAMETERS (cm)

102	87	64	36	46	100
62	17.5	49	92×28×89	N.W.12.5kg	G.W.14.5kg

Aluminum Wheelchair Series



AU869LJP-43

- › Aluminum chair frame
- › Fixed armrest
- › Fixed footrest
- › Solid castor
- › PU rear wheel
- › With united brake

MAIN PARAMETERS (cm)

96	95.5	58	42	46	100
56	15	46	96×33×90	N.W.11.8kg	G.W.13.8kg

Always care for your health

Foshan Wheelchair



FS863LAJPF4

- › Aluminum chair frame
- › Drop back handle
- › Fixed armrest
- › Flip up footrest
- › Solid castor
- › PU rear wheel
- › With united brake

MAIN PARAMETERS (cm)

100	91	62	44.5	46	100
51	15	49	80×30×75	N.W.12.2kg	G.W.14.2kg



FS872LAJP

- › Aluminum light weight chair frame
- › Fixed armrest
- › Fixed footrest
- › Solid castor
- › PU rear wheel
- › Drop back handle
- › United brake

MAIN PARAMETERS (cm)

90	92	61	44	46	100
51	15	49	82×28×76	N.W.11.7kg	G.W.14kg



FS863LABJPF3

- › Aluminum chair frame
- › Drop back handle
- › Fixed armrest
- › Flip up footrest
- › Solid castor
- › PU rear wheel
- › With united brake

MAIN PARAMETERS (cm)

100	91	62	44.5	46	100
30	15	49	71×28×75	N.W.11.5kg	G.W.13.4kg

Aluminum Wheelchair Series

• Foshan Wheelchair



FS907LABH

Aluminum Chair Frame, Flip-up Desk Armrest, Detachable Footrest, Solid Castor, Pneumatic Mag Rear Wheel, Drop Back Handle, With Drum Brake.



FS207LABP

Aluminum Chair Frame, Flip-up Inclined Armrest, Detachable Footrest, Solid Castor, PU Mag Rear Wheel, drop Back Handle.



FS953LB

Aluminum Chair Frame, Flip-up And Height-adjustable Inclined Armrest, Detachable Footrest, Angle-adjustable Footplate, Solid Castor, Pneumatic Mag Rear Wheel, Anti-tippers.



FS907LABJPF3

Aluminum Light Weight Chair Frame, Flip-up And Height-adjustable Armrest, Fixed Footrest, Solid Castor, PU Mag Rear Wheel, drop Back Handle, With United Brake.



FS871LBJ

Aluminum Light Weight Chair Frame, Fixed Armrest& Footrest, Solid Castor, Pneumatic Mag Rear Wheel, With United Brake.



FS805LABJP

Aluminum Light Weight Chair Frame, Fixed Armrest, Flip-Up Footrest, Solid Castor, PU Mag Rear Wheel, Drop Back Handle, With United Brake.



FS976LABJPF2-43

Aluminum Light Weight Chair Frame, Fixed Armrest, Detachable Footrest, PU Mag Rear Wheel, Solid Castor, Drop BackHandle, With United Brake.



FS976LAB

Aluminum Chair Frame, Fixed Armrest, Detachable Footrest, Solid Castor, Solid Mag Rear Wheel, Drop Back Handle.

MAIN PARAMETERS (cm)	MODEL NO.	Technical Parameters														
		W	H	D	W	H	D	W	H	D	W	H	DIM	kg	NW	GW
	FS907LABH	94	63	46	40	15	51	91	40	100	74x28x79	13.4	15.5			
	FS953LB	100	62	46	40	15	50	90.5	37	100	78x33x93	14.4	16.4			
	FS207LABP	101	62	46	38	15	50	89	39	100	78x28x76	12.9	14.9			
	FS907LABJPF3	92	66	46	40.6	17.5	49	88	39	100	80x30x72	13.5	15.5			

MAIN PARAMETERS (cm)	MODEL NO.	Technical Parameters														
		W	H	D	W	H	D	W	H	D	W	H	D	kg	NW	GW
	FS871LBJ	89	59	46	40	15	49	89	39	100	80x28x90	10.8	13.2			
	FS976LABJPF2-43	96	55	43	30	20	49	93	39	100	74x25x78	12.3	13.8			
	FS805LABJP	92	59	46	31	15	48	91	40	100	71x25x76	9.9	11.5			
	FS976LAB	96	60	46	20	20	52	101	40	100	74x26x78	10.6	12			

Aluminum Wheelchair Series



15



FS800LBJF5

- › Aluminum chair frame
- › Flip-up footrest
- › PU castor
- › PU mag rear wheel
- › Flip-up armrest
- › United brake



FS800LBJF4

- › Aluminum Chair Frame
- › Flip-up Footrest
- › Solid Castor
- › Solid Mag Rear Wheel
- › Flip-up Armrest
- › United Brake



FS800LBF3

Transit Chair, Strong And Durable Aluminym Frame, Easy To Fold, 5" Rear Castors With Pedal Lock, Flip-up Armrest



FS804LABJP

Transit Wheelchair, Aluminum Light Weight Chair Frame, 8"PU Mag Rear Wheel, 5"Solid Castor, Flip-up Armrest And Footrest, Drop Back Handle, With United Brake And Seat Belt.



FS800LBJ

Transit Wheelchair, Aluminum Light Weight Chair Frame, 6" Solid Mag Rear Wheel, 5"Solid Castor, Flip-up Armrest, With United Brake And Seat Belt.



FS800LBF2

Aluminum Chair Frame, Flip-up Footrest, Solid Castor, Solid Mag Rear Wheel, Flip-up Armrest, United Brake,

MAIN PARAMETERS (cm)	MODEL NO.	90	60	45	30	15	51	95	47	100	59x32x83	10.5	11.5
		FS800LBJF5	FS800LBJF4	FS800LBF3	FS804LABJP	FS800LBJ	FS800LBF2				DIM (cm)	NW (kg)	GW (kg)
		90	60	45	30	15	51	95	47	100	59x32x83	10.5	11.5
		FS800LBJF5	FS800LBJF4	FS800LBF3	FS804LABJP	FS800LBJ	FS800LBF2				DIM (cm)	NW (kg)	GW (kg)
		84	48	40	20	15	48	89	48	100	50.5x26x81	7.5	9
		73.5	38	25	12.5	7.6	49	86.5	34	100	42x30x61	6	8
		84	45	37	20	12.5	45	88	39.5	100	75x25x65	7.1	9.4
		68	51	38	15	12.5	45	88	33	100	48x34x66	7	9.1
		83	50	38	20	13	49	90	34	100	51.5x25.5x76.5	7	9.2

16



Aluminum Wheelchair Series



FS874LAJ-30

- › Pediatric wheelchair
- › Aluminum chair frame
- › Fixed armrest
- › Fixed footrest
- › Solid castor
- › Pneumatic rear wheel
- › Drop back handle
- › With united brake

MAIN PARAMETERS (cm)

85	86	51	40	30	75
53	15	44	78x28x68	N.W.11.4kg	G.W.14.4kg

17



FS980LQF3-30

- › Aluminum pediatric wheelchair
- › Quick release pneumatic rear wheel
- › Flip up, height adjustable armrest
- › Detachable footrest
- › Solid castor

MAIN PARAMETERS (cm)

85	88	46	32	30	75
52	15	52	74x33x91	N.W.12.5kg	G.W.15.5kg

18



AU974LA-35

- › Pediatric wheelchair
- › Aluminum chair frame
- › Fixed armrest
- › Detachable footrest
- › Drop-back handle
- › Pneumatic rear wheel
- › Solid castor
- › Anti-tippers

MAIN PARAMETERS (cm)

93	87	55	37	35	100
53	10	42	69x33x66	N.W.13.8kg	G.W.16kg

Aluminum Wheelchair Series



FS980LA-46

- › Aluminum chair frame
- › Flip-up armrest
- › Fixed footrest
- › Solid castor
- › Pneumatic rear wheel
- › Drop back handle

MAIN PARAMETERS (cm)

W	D	H	AS	AF	SH
103	94	64	40	46	100
61	15	51	94x33x75	DIM	N.W.14.6kg G.W.17.6kg

Always care for your health



FS980LA-35

- › Pediatric wheelchair
- › Aluminum chair frame
- › Flip-up armrest
- › Fixed footrest
- › Solid castor
- › Pneumatic rear wheel
- › Drop back handle

MAIN PARAMETERS (cm)

W	D	H	AS	AF	SH
93	87	55	37	35	100
53	10	42	69x33x66	DIM	N.W.13.8kg G.W.16kg

Foshan Wheelchair



FS958LBPYF8

- › Aluminum pediatric wheelchair
- › Adjustable seat size and seat height
- › Height adjustable footrest and adjustable headrest
- › Movable armrest
- › Optional accessories: food table
- › Tiltable seat
- › Reclining back
- › PU wheels
- › Push handle

MAIN PARAMETERS (cm)

W	D	H	AS	AF	SH
75	98-120	52	32	25-32	35
30	15	38-66	86x53x45	DIM	N.W.24kg G.W.27kg



FS980LQF8

- › Aluminum pediatric wheelchair
- › Fixed footrest, height adjustable footplate
- › Quick release pneumatic rear wheel
- › Adjustable seat size
- › Solid castor
- › Flip-up armrest

MAIN PARAMETERS (cm)

W	D	H	AS	AF	SH
74	87-92	50-60	33	25-35	75
52	10	50-55	76x39x90	DIM	N.W.14.4kg G.W.17.4kg

Reclining Wheelchair Series



FS958LBCGJPY

- › Pediatric wheelchair
- › Detachable armrest
- › Detachable elevating footrest
- › Solid castor
- › PU mag rear wheel
- › Tilting seat
- › Adjustable headrest
- › Safety belt

MAIN PARAMETERS (cm)

113	93-113	54	42	32-38	75
38	15	48.5	77x40x95	N.W.22.2kg	G.W.25.7kg

Always care for your health



FS258LBGPY

- › Aluminum chair frame
- › Reclining high back
- › Solid castor
- › PU mag wheel
- › Tilting seat
- › Detachable seat and back
- › Safety belt

MAIN PARAMETERS (cm)

94	103	44	31	30	75
30	15	49	80x41x69	N.W.17.3kg	G.W.20.7kg

21

22



FS212BCEG

- › Chromed steel frame
- › Reclining high back
- › Flip-up desk armrest
- › Elevating footrest
- › Solid castor
- › Pneumatic mag rear wheel
- › Detachable push bar

MAIN PARAMETERS (cm)

106	107	55	58	41	100
40	15	49	83x33x109	N.W.18.7kg	G.W.21.4kg

Foshan Wheelchair



FS985LBGJPY

- › Aluminum chair frame
- › Reclining high back
- › Adjustable and detachable head rest
- › Detachable armrest
- › Height adjustable armrest pad
- › Solid castor
- › Height-adjustable footplate
- › PU mag rear wheel
- › Safety belt
- › United brake

MAIN PARAMETERS (cm)

97	93-113	54	42	32-38	75
38	15	48.5	84x40x95	N.W.18.7kg	G.W.21.4kg

Reclining Wheelchair Series



FS958LABGJPF10

- › Aluminum chair frame
- › Detachable and height adjustable armrest
- › Detachable footrest
- › Reclining and foldable back
- › Tilting seat
- › With united brake

MAIN PARAMETERS (cm)

123	103-130	57	54	44	75
30	15	55	80×35×89	N.W.20.8kg	G.W.29.7kg

Always care for your health



AU004LBCGJP

- › Aluminum chair frame
- › Reclining high back
- › Adjustable and detachable head rest
- › Flip-up and detachable armrest
- › Height adjustable armrest pad
- › Elevating footrest
- › Solid castor
- › PU mag rear wheel
- › United brake

MAIN PARAMETERS (cm)

117	129	67	45	47	100	
61	20	52.5	90×36×98	N.W.22.3kg	G.W.24.3kg	



FS985LAF1

- › Aluminum chair frame
- › Reclining high back
- › Drop back handle
- › Adjustable and detachable headrest
- › Flip up and height-adjustable armrest
- › Detachable footrest
- › Tilting seat
- › Solid castor
- › Pneumatic rear wheel

MAIN PARAMETERS (cm)

106	128	72	46	46	100
62	20	52	85×39×77	N.W.19.2kg	G.W.21.5kg

Foshan Wheelchair



FS954LGCJPF4

- › Aluminum chair frame
- › Reclining high back
- › Adjustable and detachable head rest
- › Flip-up and detachable armrest
- › Height adjustable armrest pad
- › Elevating footrest
- › Solid castor
- › PU mag rear wheel
- › United brake

MAIN PARAMETERS (cm)

125	122	67	43	46	100
58	15	47	96×35×85	N.W.18.9kg	G.W.22.2kg



AU204L

- > Aluminum chair frame
- > Comfortable ergonomic seat
- > Detachable and height adjustable armrest
- > Detachable elevating footrest
- > Solid castor
- > Pneumatic rear wheel

MAIN PARAMETERS (cm)

116	109-125	66.5	48	41.5	100
62	20	58	96x66x107	N.W.24.8kg	G.W.28kg

25



FS204GCBJP

- > Powder coating steel frame
- > Reclining high back
- > Adjustable and detachable headrest
- > Tilting seat
- > Height adjustable armrest
- > Elevating footrest
- > TPR castor
- > United brake
- > PU mag rear wheel
- > Optional accessories: neck cushion

MAIN PARAMETERS (cm)

121	121-133	66.5	55	46	150
61	20	63	88x64x84	N.W.38.5kg	G.W.42.8kg



FS985LAJPF2-46

Aluminum Frame, Flip-up Armrest, Height Adjustable Armrest Pad, Detachable Footrest, Pu Rear Wheel, Pu Castors, Foldable Back



FS954LGC

Aluminum Chair Frame, Reclining High Back, Adjustable And Detachable Headrest, Detachable Armrest, Elevating Footrest, Solid Castor, Pneumatic Rear Wheel, Plastic Hand Rim.



FS954GCJ-46

Powder Coating Steel Frame, Reclining High Back, Detachable Armrest Elevating Footrest, Solid Castor, Solid Rear Wheel, With United Brake Optional Accessories: Neck Cushion



FS902GC-46

Powder Coating Steel Frame, Reclining High Back, Flip Up Armrest, Elevating Footrest, Solid Castor, Solid Rear Wheel, Optional Accessories: Neck Cushion

MAIN PARAMETERS (cm)	MODEL NO.														DIM	N.W. kg	G.W. kg
	FS985LAJPF2-46	99	68	46	50.8	17.5	50	119	44	100	81x41x73	19.6	22				
	FS954LGC	124	68	46	62	20	51	128	44	100	98x35x89	18.3	21				
	FS954GCJ-46	123	65	46	58	20	51	119	45	100	98x35x89	27.7	27.7				
	FS902GC-46	123	65	46	58	20	50	119	45	100	95x34x92	23.5	26.5				

26

FS129F1

- › Power standing wheelchair
- › 4 Positioning modes:
Seating / Elevating / Reclining / Standing
- › Strong body / Waist / Kneel support



27



MAIN PARAMETERS (cm)

MODEL NO.																		
FS129F1	132	65	43	35	35	56-69	133	54	22	1-6	55AH	320Wx2	5	12°	100	118x71x86	143	179



FS139

- › Power stand up wheelchair
- › Powder coating steel frame
- › Flip-up armrest with single controller for both driving & seat lifting
- › Rubber foam castor
- › Strong rubber foam drive wheel
- › Special design with comfortable safety legrest belt

MAIN PARAMETERS (cm)

99	101.5-143.5	63	51	46	75	
20	1-7	28	200Wx2	4	12°	
32	20	52	97x64x68	N.W.71kg	G.W.85kg	

28



FS129

- › Power stand up wheelchair
- › Powder coating steel frame
- › Flip-up armrest with single controller for both driving & seat lifting
- › Rubber foam castor
- › Strong rubber foam drive wheel
- › Special design with comfortable safety legrest belt

MAIN PARAMETERS (cm)

107	107	64	48	46	100	
25	1-7	40AH	200Wx2	4	12°	
32	20	58	110x64x75	N.W.82kg	G.W.98.2kg	

Power Wheelchair Series



FS123

- › Powder coating steel frame
- › Adjustable armrest pad
- › Detachable armrest
- › Detachable footrest
- › Angle-adjustable backrest
- › Rubber foam castor and drive rear wheel
- › With safety lamps
- › Lead acid battery

MAIN PARAMETERS (cm)

107	97	64	47	46	100
20	1-7	28AH	200W×2	4	12°
32	20	55	78×68×69	N.W.63kg	G.W.70.8kg

Always care for your health



FS123GC

- › Powder coating steel frame
- › Adjustable armrest pad
- › Detachable armrest
- › Elevating footrest,
- › Rubber foam castor
- › Rubber foam drive rear wheel
- › Angle-adjustable backrest
- › With safety lamps
- › Adjustable headrest

MAIN PARAMETERS (cm)

123	93-126	67	47	43	125
25	1-8	40AH	200W×2	5	12°
35	23	58	78×68×69	N.W.81kg	G.W.88kg

29

30



FS123F1

- › Powder coating steel frame,
- › Adjustable armrest pad,
- › Detachable footrest,
- › Angle-adjustable backrest,
- › PU rear wheel,
- › PU castor,
- › With safety lamps,

MAIN PARAMETERS (cm)

112	95	67	47	43	125
25	1-8	40AH	200W×2	5	12°
35	23	53	78×68×69	N.W.77.3kg	G.W.84kg

Foshan Wheelchair



FS124

- › Powder coating steel frame,
- › Adjustable and detachable armrest pad,
- › Detachable armrest
- › Rubber foam castor and drive rear wheel.
- › Powered lifting seat,
- › Detachable footrest
- › Lead acid battery

MAIN PARAMETERS (cm)

104	92	62	39	46	100
25	1-7	40AH	200W×2	4	12°
32	20	52	76×64×63	N.W.82.8kg	G.W.91kg

Power Wheelchair Series



FS131GC

- › Steel chair frame,
- › Adjustable and detachable armrest,
- › Reclining high back,
- › Detachable and elevating footrest,
- › Rubber foam castor and drive rear wheel,
- › Electrical lifting seat,
- › Tilting seat.

MAIN PARAMETERS (cm)

108	106	64	37	46	91.7
25	1-7	40AH	200W×2	4	12°
32	20	59	85×65×75	N.W.91.7kg	G.W.98.5kg

Always care for your health



FS123GF3

- › Powder coating steel frame,
- › Adjustable and detachable pu armrest,
- › Detachable footrest,
- › Angle-adjustable footplate,
- › Pneumatic castor and rear wheel,
- › Adjustable seat angle,
- › With safety lamps,
- › Adjustable headrest,
- › Joystick can move to back handle,
- › Second seat for caretaker.

MAIN PARAMETERS (cm)

170	93-126	67	47	43	125
30-40	1-6	40AH	200W×2	5	12°
35	23	58	169×68×69	N.W.84.5kg	G.W.90kg



FS130LGC

- › Aluminum chair frame,
- › Comfortable vehicle seat,
- › Reclining high back,
- › Elevating footrest,
- › Detachable armrest,
- › Rubber foam castor and drive rear wheel.

MAIN PARAMETERS (cm)

112	92	61	55	46	100
20	1-6	28AH	200W×2	4	12°
32	20	56	77×61×83	N.W.69.8kg	G.W.75kg



FS141

- › Comfortable and rotatable seat,
- › Flip-up armrest,
- › Driving lamp,
- › Rubber foam castor and foam rear wheel.

MAIN PARAMETERS (cm)

123	129-139	58	56	48	100
30-40	1-10	33AH	400W	5	12°
26	22	58	124×59×80	N.W.78kg	G.W.90.5kg

Foshan Wheelchair



FS112LAF2

- Aluminum power wheelchair,
- Flip-up armrest,
- Flip-up footrest,
- Foldable back,
- Lithium battery,
- PU castor,
- PU rear wheel

MAIN PARAMETERS (cm)

104	91	63	49.5	46	75	
17	1-6	15.6AH	200W×2	4	12°	
32	20	49.5	80×40×77	N.W.26.4kg	G.W.30.9kg	



FS110LABF3

Aluminum Frame, Drop Back Handle, Flip-up Armrest, Detachable Footrest, Rubber Foam Castor And Rear Wheel, Lithium Battery On Both Armrest Board.



FS110LAEPF2

Aluminum Chair Frame, Drop Back Handle, Flip-up Armrest, Detachable Footrest, Rubber Foam Castor And Rear Wheel, Lithium Battery On Both Armrest Board, Foldable Push Bar



FS112LAF4

- Aluminum power wheelchair,
- Flip-up armrest,
- Flip-up footrest,
- Foldable back,
- Lithium battery,
- PU castor,
- PU rear wheel

MAIN PARAMETERS (cm)

101	92.5	62	44	46	75	
17	1-6	15.6AH	250W×2	4	12°	
41	20	49.5	82×46×77	N.W.28.5kg	G.W.33kg	



FS111AF1

Powder Coating Steel Frame, Drop Back Handle, Flip-up Pu Armrest, Detachable Footrest, PU Castor And Rear Wheel, Foldable Chair With Lead Acid Battery On Both Side.



FS110AF2

Powder Coating Steel Frame, Drop Back Handle, Flip-up Pu Armrest, Detachable Footrest, Rubber Foam Castor And Drive Rear Wheel, Foldable Chair Frame.

MODEL NO.	MAIN PARAMETERS (cm)												MAIN PARAMETERS (kg)	DIM	NW G.W		
FS110LABF3	104	63	46	32	20	52	97.5	45	20	1-6	23AH	150W×2	4	12°	100	77×38×78	33 38.5
FS110LAEPF2	104	63	46	32	20	52	96.5	45	20	1-7	23AH	200W×2	4	12°	100	83×38×73	40 46.6
FS111AF1	107	61	46	41	20	52	94.5	43	20	1-6	20AH	250W×2	4	12°	100	86×48×78	48 55
FS110AF2	105	62	46	32	20	50	93	44	20	1-7	28AH	200W×2	4	12°	100	86×40×76	60 65



Power Wheelchair Series



FS101A

- › Powder coating steel frame,
- › Detachable armrest,
- › Detachable footrest,
- › Rubber foam castor,
- › Rubber foam rear wheel,
- › Drop back handle,
- › Foldable frame,
- › Anti tippers,

MAIN PARAMETERS (cm)

106	97	70	47	46	100
20	1-7	28AH	200W×2	4	8°
53	20	53	86×46×78	N.W.60kg	G.W.67.8kg

Always care for your health



FS111APF2

- › Powder coating steel frame,
- › Drop back handle,
- › Flip-up armrest,
- › Detachable footrest,
- › PU rear wheel,
- › Foldable chair frame with lead acid battey on both sides.

MAIN PARAMETERS (cm)

108	94	65	43	46	100
20	1-6	21AH	150W×2	4	8°
54	20	51	86×48×78	N.W.50kg	G.W.54kg



FS108LAEP

- › Aluminum chair frame,
- › Flip-up inclined armrest,
- › Detachable footrest,
- › PU castor,
- › PU rear wheel,
- › Lithium battery,
- › Foldable push bar.

MAIN PARAMETERS (cm)

107	94	68	41	46	75
20	1-7	23AH	150W×2	4	8°
60	20	49	86×41×73	N.W.32.5kg	G.W.36kg

Foshan Wheelchair



FS101LAEPF1

- › Aluminum foldable frame,
- › Drop back handle,
- › Flip-up armrest,
- › Detachable footrest,
- › PU castor,
- › PU rear wheel with hand rim,
- › With additional cushion,
- › Lithium battery on side board,
- › Foldable push bar,

MAIN PARAMETERS (cm)

111	96	68	44	46	100
20	1-6	23AH	150W×2	4	12°
60	20	53	93×42×73	N.W.33.4kg	G.W.39kg

Power Wheelchair Series



FS102LAEF2

- › Aluminum frame,
- › Drop back handle,
- › Flip-up pu armrest,
- › Flip-up footrest,
- › Rubber foam castor and drive rear wheel,
- › Foldable chair frame.

MAIN PARAMETERS (cm)

85.5	97	63	41	46	75
17	1-6	15.6	200W×2	4	12°
32	20	52	73×45×84	N.W.25kg	G.W.29.3kg

Always care for your health



FS112LAF4-35

- › Power wheelchair for children,
- › Flip-up armrest,
- › Height-adjustable footrest,
- › Foldable back

MAIN PARAMETERS (cm)

84	88.5	52.5	38	35	75
8	1-6	7.8AH	250W×2	4	8°
32	15	44	82×46×75	N.W.24.1kg	G.W.28.6kg



FS121L

- › Aluminum chair frame
- › Fixed armrest
- › Flip-up footrest
- › Lithium battery

MAIN PARAMETERS (cm)

103	90.5	60	42	42	75
15	1-6	15.6	250W×2	4	12°
32	20	48	61×41×94	N.W.21kg	G.W.25kg



FS121LF1

- › Power pediatric wheelchair,
- › Aluminum chair frame,
- › Adjustable pu armrest,
- › Adjustable seat width,
- › Reclining high back,
- › Lithium battery,
- › Rubber foam drive wheel.

MAIN PARAMETERS (cm)

96	92	56	40	38	75
20	1-7	23AH	200W×2	4	12°
32	20	52	86×58×75	N.W.40kg	G.W.47kg

Foshan Wheelchair



Leisure & Sports Wheelchair Series



FS771LQ-32

- › Dedicated wheelchair for disabled rugby matches
- › Aluminum frame
- › Fixed armrest
- › Fixed footrest
- › Solid front wheel
- › Quick release light weight pneumatic rear wheel
- › Additional sponge seat and back cushion
- › Anti tipper
- › With seat belt , leg belt and foot belt

FS770LQ-32

- › Dedicated wheelchair for disabled rugby matches
- › Aluminum frame
- › Fixed armrest
- › Fixed footrest
- › Solid front wheel
- › Quick release light weight pneumatic rear wheel
- › Additional sponge seat and back cushion
- › Anti tipper
- › With seat belt , leg belt and foot belt



FS779LQ-36 (Basketball Guard)

Dedicated Wheelchair For Disabled Basketball Matches, Aluminum Frame, Fixed Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt, Leg Belt And Foot Belt.



FS778LQ-36 (Basketball Forward)

Dedicated Wheelchair For Disabled Basketball Matches, Aluminum Frame, Fixed Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt, Leg Belt And Foot Belt.



FS777LQ-36 (Basketball Center)

Dedicated Wheelchair For Disabled Basketball Matches, Aluminum Frame, Fixed Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt , Leg Belt And Foot Belt.



AU778LQ-36

Dedicated Wheelchair For Disabled Basketball Matches, Aluminum Frame, Height-adjustable Back, Fixed Armrest, Fixed Footrest, Height-adjustable And Level-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt, Leg Belt And Foot Belt.



FS710L-30

Dedicated Wheelchair For Sports Events Such As Racing And Road Marathons In Disabled Venues, Aluminum Frame, 18" PU Front Wheel With Aluminum Brake, 28" Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Seat Belt Optional Seat Width:28cm/30cm/32cm/34cm/36cm.

41



FS756LQ-36

Dedicated Wheelchair For Disabled Ping Pong Matches, Aluminum Frame, Fixed Armrest, Fixed Footrest, Height-adjustable And Level-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt .



FS766LQ-36

Dedicated Wheelchair For Disabled Fencing Matches, aluminum Frame, Fixed Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Seat Belt .

42



FS755LQ-36

Dedicated Wheelchair For Disabled Dancing Matches, Aluminum Frame, Fixed Armrest, Fixed Footrest, Height-adjustable And Level-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt , Leg Belt And Foot Belt.



FS786LQ-36

Dedicated Wheelchair For Disabled Badminton Matches, Aluminum Frame, Fixed Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt .



FS785LQ-36

Dedicated Wheelchair For Disabled Tennis Matches, Aluminum Frame, Fixed Armrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, Anti Tipper, Seat Belt .



AU736LQ-36

Aluminum Foldable Frame, Foldable Back, Detachable Armrest, Fixed Footrest, foldable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS730LQF10-36

Aluminum Frame, Foldable Back, Detachable Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion, shock Springs Under The Seat.



FS737LQ-36

Aluminum Frame, Foldable Back, Detachable Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS735LQ-36

Aluminum Frame, Flip-up Armrest, Fixed Footrest, Height-adjustable Footplate, Anti-skidding Solid Front Wheel, Quick Release Anti-skidding Pneumatic Rear Wheel.



FS734LQ-36

Aluminum Frame, Foldable Back, Detachable Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS730LQ-36

Aluminum Frame, Foldable Back, Detachable Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS731LQ-36

Aluminum Frame, Foldable Back, Detachable Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS732LQ-36

Aluminum Frame, Foldable Back, Fixed Armrest, Fixed Footrest, Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.

MAIN PARAMETERS (cm)	MODEL NO.	Dimensions and Weight											
		Width (cm)	Length (cm)	Height (cm)	Seat Width (cm)	Seat Depth (cm)	Armrest Height (cm)	Footrest Height (cm)	Wheelchair Weight (kg)	Front Wheel Diameter (cm)	Rear Wheel Diameter (cm)		
	AU736LQ-36	85	52	36	60	10	52	82	32	100	86x38x67	13.4	15
	FS737LQ-36	86	64	36	60	10	47	79	35	100	84x61x70	10.5	14
	FS730LQF10-36	76	51	36	61	7.5	54	75	24	100	78x52x63	11.5	15
	FS735LQ-36	81	55	36	60	10	49	84	31	100	83x34x84	10.3	13

MAIN PARAMETERS (cm)	MODEL NO.	Dimensions and Weight											
		Width (cm)	Length (cm)	Height (cm)	Seat Width (cm)	Seat Depth (cm)	Armrest Height (cm)	Footrest Height (cm)	Wheelchair Weight (kg)	Front Wheel Diameter (cm)	Rear Wheel Diameter (cm)		
	FS734LQ-36	90	59	36	60	10	51	78	32	100	79x60x59	12.6	15.9
	FS731LQ-36	86	65	36	60	10	52	86	34	100	76x62x63	10.2	12
	FS730LQ-36	85	67	36	60	10	48	78	33	100	76x62x50	9.5	11.8
	FS732LQ-36	89	58	36	60	10	50	76	33	100	76x62x50	9.3	12



FS727LQF1-36

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS722LQ-36

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable And Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Rear Wheels Can Be Adjusted Forward And Backward,additional Sponge Seat And Back Cushion.



FS727LQ-36

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable And Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS720LQ-36

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable And Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS723LQF15-38

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS723LQF1-36

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable And Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion,anti Tipper, Luggage Net Under Seat.



FS723LQ-36

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable And Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.



FS721LQ-36

Aluminum Foldable Frame, Fixed Armrest, Fixed Footrest, Foldable And Height-adjustable Footplate, Solid Front Wheel, Quick Release Light Weight Pneumatic Rear Wheel, Additional Sponge Seat And Back Cushion.

MAIN PARAMETERS (cm)	MODEL NO.	宽度 (cm)	深度 (cm)	高度 (cm)	座宽 (cm)	座深 (cm)	扶手高度 (cm)	座面高度 (cm)	座宽 (cm)	座深 (cm)	扶手高度 (cm)	座面高度 (cm)	重量 (kg)	尺寸 (cm)	净重 (kg)	毛重 (kg)
	FS727LQF1-36	81	55	36	60	10	49	84	31	100	83x34x84	10.3	13			
	FS722LQ-36	84	57	36	60	10	48	72	28	100	89x36x76	11.6	13.6			
	FS723LQF15-38	85	54.5	38	60	10	50	80	27	100	85x34x83	12.9	15.1			
	FS723LQF1-36	84	60	38	60	10	49	77	30	100	87x36x79	12.8	15.2			

MAIN PARAMETERS (cm)	MODEL NO.	宽度 (cm)	深度 (cm)	高度 (cm)	座宽 (cm)	座深 (cm)	扶手高度 (cm)	座面高度 (cm)	座宽 (cm)	座深 (cm)	扶手高度 (cm)	座面高度 (cm)	重量 (kg)	尺寸 (cm)	净重 (kg)	毛重 (kg)
	FS727LQ-36	83	52	36	60	10	48	76	30	100	87x36x80	10.6	13.6			
	FS720LQ-36	83	59	36	60	10	48	72	29	100	87x36x73	11.3	13			
	FS723LQ-36	83	56	36	60	10	48	76	32	100	87x36x79	8	11			
	FS721LQ-36	83	58	36	60	10	48	72	29	100	87x36x73	11.6	13.6			



Steel Wheelchair Series

47



FS951B-56

- > Powder coating steel frame
- > 22" seat width
- > Flip-up armrest
- > Detachable footrest
- > PU castor
- > Pneumatic mag rear wheel

MAIN PARAMETERS (cm)

110	93	73	41	56	125
62	20	54	83x33x93	N.W.21kg	G.W.24.1kg



FS210ABE-61

- > Powder coating steel frame,
- > Detachable desk armrest,
- > Drop back handle,
- > Detachable footrest,
- > 24" seat width,
- > Detachable push bar,
- > PU castor,
- > Pneumatic mag rear wheel.

MAIN PARAMETERS (cm)

103	92	78	41	61	125
30	20	50	80x32x83	N.W.18.5kg	G.W.20.6kg

48



FS209AE-61

- > Powder coating steel frame
- > Detachable desk armrest
- > Detachable footrest
- > 24" seat width
- > Detachable push bar
- > PU castor
- > Pneumatic rear wheel
- > Drop back handle

MAIN PARAMETERS (cm)

111	92	81	41	61	125
62	20	52	86x33x83	N.W.22kg	G.W.24.5kg

Steel Wheelchair Series



FS975-51

- > Chromed steel frame
- > 20" seat width
- > Detachable armrest
- > Detachable footrest
- > Solid castor
- > Pneumatic rear wheel

MAIN PARAMETERS (cm)

105	91	73	40	51	125
62	20	50	80x28x93	DIM	N.W.18.7kg G.W.21.3kg

Always care for your health



FS908AQ

- > Chromed steel frame
- > Flip-up desk armrest
- > Detachable footrest
- > Quick release pneumatic rear wheel
- > Drop back handle
- > Solid castor

MAIN PARAMETERS (cm)

106	91	64	42	46	100	
61.5	15	51	80x28x78	DIM	N.W.15.8kg G.W.18.6kg	



FS211BF1

- > Powder coating steel frame
- > Detachable armrest
- > Detachable & length-adjustable footrest
- > Solid mag rear wheel
- > Solid castor

MAIN PARAMETERS (cm)

107	91.5	67	40	46	100	
62	20	49	82x31x93	DIM	N.W.19.4kg G.W.21.9kg	



FS909B

- > Powder coating steel frame
- > Flip-up armrest
- > Detachable footrest
- > Solid castor
- > Solid mag rear wheel

MAIN PARAMETERS (cm)

107	93	62	40	46	100	
62	20	49	82x28x95	DIM	N.W.18.2kg G.W.20.5kg	

Foshan Wheelchair



FS909

- › Powder coating steel frame
- › Flip-up desk armrest
- › Detachable footrest
- › Solid castor
- › Solid rear wheel

MAIN PARAMETERS (cm)

105	92	63	39	46	100
58	20	50	82×28×95	N.W.16.9kg	G.W.19.4kg



FS902C

Chromed Steel Frame, Detachable Desk Armrest, Elevating Footrest, TPR Castor, Solid Rear Wheel.



FS903

Powder Coating Steel Frame, Detachable Armrest, Detachable Footrest, Solid Castor, Solid Rear Wheel.



FS909C

- › Chromed steel frame
- › Flip-up armrest
- › Detachable-elevating footrest
- › Solid rear wheel
- › Solid castor

MAIN PARAMETERS (cm)

109	92	63	39	46	100
58	20	50	82×28×93	N.W.18.4kg	G.W.20.9kg

FS901B

Powder Coating Steel Frame, Detachable Armrest, Detachable Footrest, PU Castor, Pneumatic Mag Rear Wheel.



FS901

Chromed Steel Frame, Detachable Armrest, Detachable Footrest, PU Castor, Pneumatic Rear Wheel.

MAIN PARAMETERS (cm)	MODEL NO.	Dimensions (cm)							DIM	NW (kg)	GW (kg)	
FS902C	107	63	46	58	20	50	87	40.5	100	80×28×89	19.1	22
FS903	107	65	46	58	20	50	86	40	100	80×28×89	17.9	19.5
FS901B	107	64	46	61	20	50	88	40	100	80×28×89	18.5	21
FS901	107	65	46	62	20	51	88	40	100	80×28×89	18.4	20.4

Steel Wheelchair Series

—• Foshan Wheelchair



FS972

Chromed Steel Frame, Fixed Armrest, Detachable Footrest, Solid Castor, Solid Rear Wheel.



FS972B

Powder Coating Steel Frame, Fixed Armrest Detachable Footrest, Solid Castor, Solid Mag Rear Wheel.



FS873

Chromed Steel Frame, Powder Coating Cross Bar, Fixed Armrest, Fixed Footrest, Solid Castor, Solid Rear Wheel.



FS810Y

Chromed Steel Frame, Double Cross Bar, Fixed Armrest, Fixed Footrest, Solid Castor, Solid Rear Wheel, Hard Cushion Seat.

53



FS875

Powder Coating Steel Frame, Fixed Armrest, Fixed Footrest, Solid Castor, Solid Rear Wheel.



FS874

Chromed Steel Frame, Double Cross Bar, Fixed Armrest,
Fixed Footrest,solid Castor, Solid Rear Wheel.

A standard manual wheelchair with a blue padded seat and backrest, and a light-colored plastic footrest.

FS869X

Powder Coating Steel Frame, Double Cross Bar, Fixed Incline Armrest, Fixed Footrest, Solid Castor, Solid Rear Wheel.

A side-view photograph of a manual wheelchair. The chair has a dark blue padded seat and backrest. The frame is made of shiny silver-colored metal. It features two large front wheels and two smaller rear wheels. The footrests are also made of dark blue plastic. The wheelchair is positioned at an angle, facing towards the left.

FS809F7

Chromed Steel Frame, Fixed Armrest, Detachable Footrest, Detachable Backrest, Solid Castor, Solid Rear Wheel.

	MODEL NO.									DIM	N.W. KG	G.W. KG	
MAIN PARAMETERS (cm)	FS972	107	63	46	58	20	49	87	40	100	80×23×89	17	19.2
	FS972B	107	63	46	62	20	49	89	40	100	80×25×91	17.8	20
	FS875	105	65	46	58	20	50	86.5	38	100	94×28×89	15.9	18.6
	FS874	106	65	46	58	20	47	92	40	100	94×28×93	16.8	19.4

	MODEL NO.												
MAIN PARAMETERS (cm)	FS873	106	62	46	58	20	50	86	38	100	94×23×89	16.2	18.2
	FS810Y	105	66	46	58	20	50	89	38	100	94×31×94	19.7	22.2
	FS869X	100	62	46	58	20	51	87	40	100	103×25×89	15.7	19
	FS809F7	106	64	46	58	20	49	87	40	100	79×22×79	16	17.5

Steel Wheelchair Series

—• Foshan Wheelchair



FS809J-35

Powder Coating Steel Frame, Fixed Armrest, Fixed Footrest, Solid Rear Wheel, Solid Castor, With United Brake



FS809B

Chromed Steel Frame, Fixed Armrest, Fixed Footrest
Solid Castor, Solid Mag Rear Wheel.



FS909P-3

Powder Coating Steel Frame, Double Cross Bar, Flip-up Arms, Detachable Footrest, Solid Castor, PU Rear Wheel.



FS976ABJ-43

Chromed Steel Frame, Fixed Armrest, Detachable Footrest, Drop Back Handle With United Brake, Solid Wheel.



FS809J

Chromed Steel Frame, Fixed Armrest, Fixed Footrest,
Solid Castors, Solid Rear Wheels, With United Brake



FS809Y

Chromed Steel Frame, Fixed Armrest, Fixed Footrest,
Solid Castors, Solid Rear Wheels, Hard Cushion Seat



FS907BC

Chromed Steel Frame, Detachable Desk Armrest,
Elevating Footrest, Solid Castor, Pneumatic Mag Rear Wheel



FS904B

Powder Coating Steel Frame, Detachable Desk Armrest, Detachable Footrest, Solid Castor, Pneumatic Mag Rear Wheel.

Model No.									DIM	N.W	G.W
MAIN PARAMETERS (cm)	FS809J-35	96	54	35	88	15	46	88	38	75	89x23x89
	FS809B	106	65	46	62	20	49	89	38	100	94x25x91
	FS809J	102	63	46	58	20	49	87	39	100	92x21x88
	FS809Y	103	64	46	58	20	50	87	38.5	100	95x31x95

	MODEL NO.									DIM	NW KG	GW KG	
MAIN PARAMETERS (cm)	FS909P-31	80	47	31	52	15	45	88	40	75	71×26×90	13.5	15.5
	FS976ABJ-43	97	56	43	20	20	50	96	42	100	74×26×78	12.9	15.1
	FS907BC	97	58	41	40	20	50	88	40	100	73×28×89	16.1	18.3
	FS904B	94	62	46	31	20	50	88	40	100	75×25×89	14.3	16.4

Steel Wheelchair Series



FS863AJPF1

- › Powder coating frame
- › Fixed armrest
- › Flip-up footrest
- › Solid castor
- › PU rear wheel
- › Drop back handle
- › United brake

MAIN PARAMETERS (cm)

96	91	63	41	46	100
52	15	48	86x31x76	N.W.14.4kg	G.W.16.4kg

Always care for your health



FS875J

- › Powder coating steel frame
- › Fixed armrest
- › Fixed footrest
- › Solid castors
- › Solid rear wheels
- › With united brake

MAIN PARAMETERS (cm)

104	87	66	40.5	46	100
58	20	48	94x25x89	N.W.15.8kg	G.W.17.8kg



FS863ABJ

- › Powder coating steel frame
- › Double cross bar
- › Drop back handle
- › United brake
- › Pneumatic mag rear wheel
- › Solid castor
- › Fixed armrest
- › Flip-up footrest

MAIN PARAMETERS (cm)

96	91	63	41	46	100
30	15	48	76x30x76	N.W.13.1kg	G.W.14.4kg

Foshan Wheelchair



FS809

- › Chromed steel frame
- › Fixed armrest
- › Fixed footrest
- › Solid castor
- › Solid rear wheel

MAIN PARAMETERS (cm)

106	87	64	38	46	100
58	20	49	92x21x88	N.W.16.7kg	G.W.19kg



Commode Series

59



FS609GCU

- › Chromed plating steel frame
- › U shape commode seat
- › Reclining back
- › Detachable armrest
- › Elevating footrest
- › Solid rear wheels , solid castors

MAIN PARAMETERS (cm)

115	129	64	79	44	100
58	20	53	107×32×80	N.W.25.6kg	G.W.29.7kg



FS654LGCJPU

- › Aluminum chair frame,
- › Detachable armrest,
- › Detachable elevating footrest,
- › U shape commode seat,
- › PU rear wheel, solid castor.

MAIN PARAMETERS (cm)

128	122	66	42	46	100
61	20	52	97×41×86	N.W.21.1kg	G.W.23.5kg

60



FS6952L

- › Aluminum chair frame
- › Adjustable headrest
- › Flip-up armrest
- › Detachable footrest
- › Tilting seat
- › Solid castors with brake

MAIN PARAMETERS (cm)

90	105-132	88	45	46	100
15	15	50	86×31×76	N.W.14.4kg	G.W.16.4kg

Commode Series



FS683L

- › Aluminum chair frame,
- › PU commode seat,
- › Flip-up armrest,
- › Detachable footrest,
- › Solid castor,
- › Pneumatic rear wheel.

MAIN PARAMETERS (cm)

106	58	61	40	46	100
58	15	48	80×35×71	N.W.15kg	G.W.17.4kg



FS6971S

- › Stainless steel frame,
- › Foldable chair frame,
- › Flip-up armrest,
- › Detachable footrest,
- › Solid castor with brake,
- › Solid mag rear wheel,

MAIN PARAMETERS (cm)

104	94	63	39	41	100
60	12.5	52	95×25×100	N.W.16.8kg	G.W.19kg



FS809LJPU

- › Chromed plating steel frame
- › Fixed armrest
- › Fixed footrest
- › U shape commode seat
- › Solid castor
- › PU rear wheel
- › With united brake

MAIN PARAMETERS (cm)

103	88	66	34.5	46	100
61	20	55	94×29×89	N.W.16.4kg	G.W.19kg



FS681U

- › Foldable chromed steel frame,
- › U shape commode seat,
- › Flip-up armrest and detachable footrest,
- › Solid castor and solid rear wheel.

MAIN PARAMETERS (cm)

102	90	60	37	45	100
58	15	53	76×29×91	N.W.19.6kg	G.W.21.4kg

Foshan Wheelchair

Commode Series



FS609U

- › Foldable chromed steel frame,
- › U shape commode seat,
- › Fixed armrest and footrest,
- › Solid castor and rear wheel.

MAIN PARAMETERS (cm)

Width	Length	Height	Seat Width	Seat Depth	Total Height
105	87	64	37	44	100
58	20	50	92×29×89	DIM	N.W.19kg G.W.21.3kg



FS609LBJPUF9-44

- › Aluminum chair frame,
- › Fixed armrest,
- › Fixed footrest,
- › PU rear wheel,
- › Solid castor,
- › With united brake

MAIN PARAMETERS (cm)

Width	Length	Height	Seat Width	Seat Depth	Total Height
89	91	58.5	37	44	100
38	15	50	88×29×93	DIM	N.W.13kg G.W.15.5kg



FS609LBJPUF10-44

- › Aluminum chair frame
- › Fixed armrest with handle
- › Fixed footrest
- › PU rear wheel
- › Solid castor
- › With united brake

MAIN PARAMETERS (cm)

Width	Length	Height	Seat Width	Seat Depth	Total Height
89	91	58.5	37	44	100
38	15	50	88×29×93	DIM	N.W.14kg G.W.16.5kg



FS609LU

- › Foldable aluminum chair frame,
- › U shape commode seat,
- › Fixed armrest and footrest,
- › Solid castor,
- › Pneumatic rear wheel.

MAIN PARAMETERS (cm)

Width	Length	Height	Seat Width	Seat Depth	Total Height
102	90	68.5	39	46	100
58	15	51	94×29×92	DIM	N.W.15.2kg G.W.18.9kg

Foshan Wheelchair

Commode Series

• Foshan Wheelchair



FS698L-A

Aluminum Commode Chair, Detachable Frame And Footrest



FS609LUF4

Aluminum Chair Frame, Adjustable Seat Height, Adjustable Backrest Height, Flip Up Armrest, Detachable Footrest, Solid Front Castor With Brake, Solid Rear Wheel,



FS697S

Stainless Steel Chair Frame, Pu Commode Seat, Flip Down Armrest, Detachable Footrest, Pedal Locking Castors.



FS697S-1

Stainless Steel Frame, Flip Up Armrest, Pedal Locking Castors, Plastic Commode Seat.



FS7964L-W5

Aluminum Chair Frame, Flip-up Armrest, Flip Up Footrest, Solid Castor, Solid Rear Wheel With Brake.



FS698L

Aluminum Chair Frame, Flip Up Footrest, Lockable Castors.



FS6951S

Stainless Steel Frame, Foldable Footrest, Lockable Castors, Flip-up Armrest.



FS695S

Stainless Steel Chair Frame, Flip Down Armrest, Pedal Locking Rear Castors, Pu Commode Seat With Cover.

MAIN PARAMETERS (cm)	MODEL NO.	DIM												NW KG	GW KG
		W	D	H	W	D	H	W	D	H	W	D	H		
	FS698L-A	97	60	49	20	15	50	99	49	100	64×33×52	10.36	13		
	FS609LUF4	94	60	42	17.5	10	47	92	27-38	100	67×60×66	12.4	14.8		
	FS7964L-W5	69	57	48	12.5	12.5	53	85	35	100	107×32×57	71	8.6		
	FS698L	84	56	45	12.5	12.5	46	89	38	100	49×24×67	8.2	10		

MAIN PARAMETERS (cm)	MODEL NO.	DIM												NW KG	GW KG
		W	D	H	W	D	H	W	D	H	W	D	H		
	FS697S	96	54	44	12.5	12.5	51	96.5	39	100	54×52.5×56	12.3	13.8		
	FS697S-1	86	55	46	12.5	12.5	54	96	38	100	56×56×58	11	12.5		
	FS6951S	91	58	41	50	20	51	93	36	100	61×55×55	12.2	14		
	FS695S	79	54	44	12.5	12.5	50	97	40	100	61×55×55	12.2	14.4		

Commode Series

• Foshan Wheelchair



FS699L

Aluminum Chair Frame, Fixed Armrest, Foldable Footrest, Solid Rear Castor With Lock, Plastic Commode Seat.



FS6992L

Aluminum Chair Frame, PU Commode Seat, Pedal Locking Rear Castors, Flip Down Armrest, Foldable Footrest.



FS691

Powder Coating Steel Frame, Detachable Armrest And Foot Rest, Solid Rear Castor With Lock. Plastic Commode Seat.



FS692

Powder Coating Steel Frame, Flip Down Armrest, Detachable Footrest, Solid Rear Castor With Lock, Plastic Commode Seat.



FS691S

Stainless Steel Chair Frame, Detachable Armrest And Footrest, Rear Wheel With Lock, Plastic Commode Seat.



FS694S

Stainless Steel Frame, PU Armrest, Solid Castor, Rear Wheel With Lock, Plastic Commode Seat, Length Adjustable Footrest.



FS692L

Aluminum Chair Frame, Flip-up Armrest, Detachable Footrest, Rear Wheel With Lock, PU Commode Seat.



FS693

Powder Coating Steel Frame, Fixed Armrest, Foldable Footrest, Rear Wheel With Lock, Plastic Commode Seat.

MAIN PARAMETERS (cm)	MODEL NO.																								
		W	D	H	W	D	H	W	D	H	W	D	H	W	D	H	W	D	H	W	D	H	W	D	H
	FS699L	67	55	45	7.5	7.5	49	92	43	100	56x56x72	8	9.6												
	FS6992L	72	55	44	12.5	12.5	55	97	40	100	98x53x54	8.4	9.7												
	FS691S	96	55	44	12.5	12.5	51	88	27	100	54x53x56	10.5	12.7												
	FS694S	83	54	44	12.5	12.5	51	84	29	100	54x53x56	9.4	11.6												

MAIN PARAMETERS (cm)	MODEL NO.																								
		W	D	H	W	D	H	W	D	H	W	D	H	W	D	H	W	D	H	W	D	H	W	D	H
	FS691	96	55	44	12.5	12.5	51	81	38	100	54x53x53	12	14												
	FS692	93	55	45	12.5	12.5	54	95	37	100	56x53x57	10.8	12.8												
	FS692L	92	55	45	12.5	12.5	54	92	41	100	56x53x57	9	10.9												
	FS693	92	54	44	12.5	12.5	47	82	38	100	50x26x66	10.5	11.7												

Commode Series

Foshan Wheelchair



FS893



FS8943



FS8964



FS8133

69



FS8997



FS893L



FS8962



FS899U

70

	MODEL NO.	寸	寸	寸	寸	寸	寸	寸	PCS/CN	G.W.(kg)	
MAIN PARAMETERS(cm)							(cm)	N.WT(kg)			
	FS893	56	44	47	94	43	38	50x19x75	7.2	1	8.5
	FS8943	60	41	47	87	41	18	46x18x56	5.3	1	6.3
	FS8997	47	46	43-53	81-91	43	42	60x52x93	15.6	2	17.4
	FS893L	56	44	47	94	43	38	50x19x75	5.3	1	6.4

	MODEL NO.	寸	寸	寸	寸	寸	寸	寸	PCS/CN	G.W.(kg)	
MAIN PARAMETERS(cm)							(cm)	N.WT(kg)			
	FS8964	52	44	45-55	78-87	43	42	52x23x79	6.2	1	7.7
	FS8133	56	44	47	85	43	38	50x19x75	7.2	1	8.5
	FS8962	51	41	44-60	76-92	42	33	58x23x78.5	6.8	1	8
	FS899U	48	43	44	78	40	35	78x44x54	4.2	4	18.8

Commode Series

Foshan Wheelchair



FS7992L



FS894L



FS8941



FS696



FS894



FS696L



FS813(S)



FS813



FS897



FS897A



FS898

	MODEL NO.	W	D	H	WD	HG	HP	PC	G.W.
MAIN PARAMETERS (cm)					(cm)	N.WT(kg)	PCS/CN		
	FS7992L	51	44	46-56	86-96	47	35	62×52×88	10.1
	FS894L	48	45	45-55	78-88	43	33	50×23×76	3.5
	FS813(S)	42	35	29-36.5	51-58.5	37	23	46×38×30	4.9
	FS813	48	48	37-47	62-72	44	33	56×40×36	6

MAIN PARAMETERS (cm)	MODEL NO.	W	D	H	WD	HG	HP	PC	G.W. (kg)
	FS8941	46	41	40-55	61-76	40	23	50×25.5×65	6.4
	FS696	55	44	7.5	7.5	45-55	78-88	50×21×76	7
	FS894	48	45	44-45	78-88	43	33	50×21×76	6
	FS696L	55	45	45-55	78-88	39	33	51×22×82	5.2
	FS897	41	38	45	40	—	—	62×30×42	9.2
	FS897A	41	38	45	40	—	—	63×25×43	5.1
	FS898	41	38	40	40	—	—	75×30×40	9.5

Commode Series

• Foshan Wheelchair



FS8951L



FS811



FS8832U



FS8992L



FS8101



FS8945



FS899-A



FS8992



FS895L



FS813L



FS8995L



FS890

MAIN PARAMETERS (cm)	MODEL NO.	寸	寸	寸	寸	寸	寸	寸	PCS CN	G.W. (kg)
									N.WT(kg)	
FS8951L	43	40	40-50	62-72	41	—	45×18×63	2.6	1	3.6
FS811	50	49	50-60	80-90	40	33	71×52×56	10	2	11.8
FS8101	46	44	41-56	60-75	37	18	46×22×70	11.5	2	12.5
FS8945	69	45.5	31-40.5	56-64.5	45	28	80×49×51	12.6	2	14.6
FS895L	51	46	42-52	70-81	43	29	60×34×52	9.2	2	10.6
FS813L	56	49	38-48	64-74	44	25	58×58×37	4.2	2	5.7

MAIN PARAMETERS (cm)	MODEL NO.	寸	寸	寸	寸	寸	寸	寸	PCS CN	G.W. (kg)
									N.WT(kg)	
FS8832U	52	45	46	85	44	38	86×16×52	8.8	1	9.8
FS8992L	60	48	45	82	41	40	88×30×60	7.6	2	8.9
FS899-A	48	43	44	78	40	35	79×43×55	16.7	4	18.7
FS8992	60	48	45	82	41	40	88×30×60	12.8	2	14.2
FS8995L	59	49	46-51	85-90	44	40	102×30×62	9.8	2	11.5
FS890	45	43	42	77	40	28	63×32×44	8	2	8.9



Bath Bench Series



FS7964L

- › Aluminum frame,
- › Flip-up armrest,
- › Flip up footrest,
- › U shape shower seat

FS7962L

- › Foldable aluminum frame,
- › Plastic seat,
- › Lockable front castor.



FS7992LU



FS790L-1



FS7985LA



FS7921LKD

MAIN PARAMETERS (cm)	MODEL NO.	Product Dimensions & Weight										PCS/CN	G.W. (kg)
		W (cm)	H (cm)	D (cm)	L (cm)	T (cm)	V (cm)	W (cm)	H (cm)	D (cm)	N.WT(kg)		
	FS7964L	63	48	51	85	38	36	104×24×52	4.7	1	5.9		
	FS7962L	75	44	48	88	39	37	67×29×53	5.9	1	7		
	FS7992LU	42	63	45-55	83-93	42	38	74×51×91.5	11	2	13.5		
	FS790L-1	50-61	56.5	36-45.5	73-82	34	39.5	57×18×87	6.67	1	8.3		
	FS7985LA	51	53	41-54	72-85	37-40	35	56×52×45	7.6	2	8.2		
	FS7921LKD	46	54	53-63	86-96	29	35	69×23×43	2.1	1	3.2		

Bath Bench Series

• Foshan Wheelchair



FS7982L



FS7991L



FS7971L



FS7963L



FS7981L



FS7986L



FS7972L



FS7985L



FS798L



FS796L



FS7987L



FS797L



FS7973L



FS790

	MODEL NO.	W	D	H	W ₁	D ₁	G	SIZE (cm)	N.WT(kg)	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)	FS7982L	42	50	35-45	74-84	30	40	68×19×39	6	2	7
	FS7991L	51	64	46-56	84-94	42	38	70×45×42	13	2	15
	FS7981L	42	50	42-52	75-85	30	40	51×47×48	8	1	7
	FS7986L	37.5	38.7	30-40	60-70	27	31.2	53×53×40	4.5	2	5.6
	FS798L	42	50	35-45	78-84	30	40	68×19×39	5.9	2	7.1
	FS796L	51	41	34-44	73-83	37	41	53×50×25	6.6	2	7.6
	FS7987L	42	50	35-45	78-84	30	40	68×19×39	5.9	2	7.1

	MODEL NO.	W	D	H	W ₁	D ₁	G	SIZE (cm)	N.WT(kg)	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)	FS7971L	42	50	42-52	45-55	30	—	50×43×45	2.6	1	3.7
	FS7963L	43	50	44-47.5	47-50.5	30	—	50×10×39	2.4	1	2.7
	FS7972L	34	32.5	38-56	38-56	32.5	—	35×35×39	3.1	2	3.9
	FS7985L	51	53	41-54	72-85	37-40	35	56×52×45	7.6	2	8.2
	FS797L	42	50	35-45	39-49	30	—	60×17×36.5	3.6	2	4.4
	FS7973L	42	50	35-45	39-49	30	—	60×17×36.5	3.6	2	4.4
	FS790	39	51	47	52	30	—	74×19×52	6.4	2	7.2

Bath Bench Series

• Foshan Wheelchair



FS7961L



FS7984LA



FS793S



FS791S



FS569S



FS568S



FS799L



FS7954



FS7931



FS794L-1



FS794L



FS7933S

MODEL NO.								N.Wt.(kg)	PCS/CN	G.W.(kg)	
MAIN PARAMETERS (cm)	FS7961L	51	41	34-44	73-83	37	41	59x49x75	7.4	2	8.8
	FS7984LA	45	42	43-53	80-90	41	37	50x46x45	3.5	1	4.7
	FS569S	51	41	14-29	29	23	—	50x42x46	6.4	2	7.6
	FS568S	27	41	14	23	41	—	45x27x19	3.2	2	3.7
	FS7931	31	66	9.5	13.5	31	—	32x7x71.5	2.8	1	3.3
	FS794L-1	23	41-46	—	18	23	—	74x25x41	9	5	10.3

MODEL NO.								N.WT(kg)	PCS/CN	G.W.(kg)	
MAIN PARAMETERS (cm)	FS793S	50	51	—	52	37	37	77×56×38.5	8.8	2	10.6
	FS791S	41	51	39-49	76-86	37	37	70×56×41.5	8.5	2	10.6
	FS799L	45	72	43-53	80-90	41	37	77×22×45	9.3	2	10.2
	FS7954	42	38	40-50	40-50	38	—	57.5×39×25	8.1	2	9
	FS794L	23	41	—	18	23	—	83×25×40	8.1	5	9.3
	FS7933S	36	84-130	—	41	28	33	74×42×27	7.2	2	7.6



Walking Aids Series



FS201

Length: 77.5cm Weight: 20kg Max Load: 200kg
Adjustable-grip angle
Armrest width adjusted from 53cm-76cm
Armrest height adjusted from 88cm-131cm
Leg width adjusted from 70cm-85cm
Strong supported by gas spring
Rear wheel with brake
Optional accessories: oxygen holder, basket, iv pole.

FS203

Length: 65cm Weight: 18kg Maxload: 120kg
Adjustable-grip angle
Armrest width adjusted from 53cm-76cm
Armrest height adjusted from 71cm-110cm
Leg width adjusted from 64cm-70cm
Strong supported by gas spring
Rear wheel with brake
Optional accessories: oxygen holder, basket, iv pole.



FS965LHF9



FS965LHF10



FS965LHF5



FS965LHF6

MAIN PARAMETERS (cm)	MODEL NO.	W	D	H	L	(cm)	N.WT(kg)	PCS/CN	G.W.(kg)
	FS965LHF9	65	62	50	81-95	66x25x83	7.3	1	9
	FS965LHF10	90	62	50	81-95	66x25x83	8.5	1	10.2
	FS965LHF5	64	61	51	86-89	50x21.5x89	8.45	1	9.9
	FS965LHF6	64	61	51	86-89	50x21.5x89	7.45	1	8.9

Walking Aids Series

• Foshan Wheelchair



FS965LHF1



FS9862LH



FS965LH-1



FS965LH

83



AU965L

Foldable Aluminum Rollater, With Hand Brake.



AU965LF1

Foldable Aluminum Rollater, With Hand Brake, And Footrest.

84



FS966LH-B



FS9865LH

	MODEL NO.	W	D	H	L	(cm)	N.WT(kg)	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)	FS965LHF1	68.5	65	51	83-87.5	69x26x86	7.8	1	8.1
	FS9862LH	65	46	55	81-93	75.5x27x84	7	1	8.7
	AU965L	68	58	56	82-102	58x33x78	7.6	1	9.4
	AU965LF1	94	58	56	82-102	58x33x78	9.3	1	11.1

	MODEL NO.	W	D	H	L	(cm)	N.WT(kg)	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)	FS965LH-1	68	46	53	78-90	62x24x84	7.4	1	9
	FS965LH	68	46	53	78-90	62x24x84	7.4	1	9
	FS966LH-B	68	46	53	78-90	61x23x86	6.6	1	8.2
	FS9865LH	63	39	43	86.5-96.5	51.5x33x79	6.17	1	8.1

Walking Aids Series

• Foshan Wheelchair



FS966LH



FS914H-1



FS9634L



FS9635L



FS9622



FS968L



FS9636L



FS961L

85

86

MAIN PARAMETERS (cm)	MODEL NO.	W	D	H	L	(cm)	N.WT.(kg)	PCS/CN	G.W.(kg)
	FS966LH	61	46	55	80-92	62×24×84	7	1	8.5
	FS914H-1	70	43	62	78-90	61×25×72	108	1	10.9
	FS9622	60	47	47	76-86	70×20×75	10.3	1	12
	FS968L	59	46	53	78-90	61×19×81	7	1	8.4

MAIN PARAMETERS (cm)	MODEL NO.	W	D	H	L	(cm)	N.WT.(kg)	PCS/CN	G.W.(kg)
	FS9634L	48	44	—	82-100	57×19×86	5.1	2	6.4
	FS9635L	48	44	—	82-100	57×19×86	5.1	2	6.4
	FS9636L	48	44	—	82-100	57×19×96	5.1	2	6.4
	FS961L	49	44	43-61	74-92	57×19×79	6.7	2	7.7

Walking Aids Series

• Foshan Wheelchair



FS9133L



FS915L



FS919L



FS967L



FS913L



FS913L-B



FS963L



FS9134L



FS917L(S/M/L)



FS9124L

MAIN
PARAMETERS
(cm)

MODEL NO.	W	D	H	L	(cm)	N.WT.(kg)	PCS/CN	G.W.(kg)
FS9133L	47	46	—	75-93	52×16×79	4.4	2	5.3
FS915L	47	47	—	75-93	53×17×79	5.1	2	6.2
FS919L	47	46	—	75-96	52×16×79	4.8	2	5.8
FS967L	48	44	—	79-96	57×18×82	5.2	2	6.5

MAIN
PARAMETERS
(cm)

MODEL NO.	W	D	H	L	(cm)	N.WT.(kg)	PCS/CN	G.W.(kg)
FS913L	47	46	—	75-93	52×16×79	4.8	2	5.7
FS913L-B	47	46	—	75-93	52×16×79	4.4	2	5.3
FS963L	47	46	—	77-95	58×18×81	5.2	2	7
FS9134L	47	46	—	76-91	52×18×51	4.8	2	5.7
FS917L(S)	43	38	—	65-75	76×68×35	5.6	4	7.6
FS917L(M)	44	38	—	75-85	86×68×35	5.7	4	7.7
FS917L(L)	45	38	—	85-95	96×68×35	5.8	4	7.8
FS9124L	35	34	—	75-82	43×14.5×76.5	3.2	2	4

Walking Aids Series

• Foshan Wheelchair



FS9631L



FS9637L



FS912L



FS9125L



FS9126L



FS9128L



FS9123L



FS969H

89

90

	MODEL NO.	寸	吋	寸	吋	寸	吋	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)									
	FS9631L	52	64	—	73.5-92	58x17x76	2.5	2	6.2
	FS9637L	48	44	—	82-100	57x19x86	5.1	2	6.4
	FS9126L	51	44	—	80-90	59x18x70	5	2	6.1
	FS9128L	55	44	—	83-93	56x18x70	7	2	8.1

	MODEL NO.	寸	吋	寸	吋	寸	吋	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)									
	FS912L	50	44	—	77-94	57x18x80	4.8	2	5.8
	FS9125L	51	44	—	77-95	62x18x82	6	2	6.8
	FS9123L	85	72	—	80-100	76x19x96	9.4	2	11.2
	FS969H	59	38	—	82-97	63x25x70	7.7	1	9.2

Walking Aids Series

• Foshan Wheelchair



FS9142L



FS914L



FS912L(XS)



FS9122L-A



FS9121L-B



FS9122L



FS9123L(S)-B



FS9145L

	MODEL NO.	W	D	H	L	(cm)	N.WT.(kg)	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)									
	FS9142L	60	43.5	55	84	55x19.5x83	5.1	1	6.2
	FS914L	50	47	—	79-89	60x34x89.5	3.5	2	8.7
	FS9121L-B	55	35	—	51-61	59x30x52	5.4	2	6.4
	FS9122L	52	59	—	45-55	66x18x65	6	2	6.9

	MODEL NO.	W	D	H	L	(cm)	N.WT.(kg)	PCS/CN	G.W.(kg)
MAIN PARAMETERS (cm)									
	FS912L(XS)	42	55	—	48.5-58.5	53x18x52	3.6	2	4.6
	FS9122L-A	52	50	—	60-70	61x25x80	8.6	2	9.7
	FS9123L(S)-B	71	—	—	54-64	68x29x59	6.6	2	7.9
	FS9145L	39	62	—	51-61	52.5x22x63	2.9	2	3.9

Walking Aids Series

• Foshan Wheelchair



FS9776



FS960



FS9600



FS9601



FS9610



FS9612



FS9611



FS9613

MODEL NO.						N.Wt(kg)	PCS/CN	G.W.(kg)	
MAIN PARAMETERS (cm)	FS9776	43	18	42-58	84-105	83x45x42	10.3	1	12.5
	FS960	58	45	—	31	51x41x34	2.7	2	5.4
	FS9610	74	45	55	110-155	74x15x85	11.9	1	13.6
	FS9612	71	44	52	109-135	77x27x97	9.4	1	11.8

MODEL NO.						N.W.(kg)	PCS/CN	G.W.(kg)	
MAIN PARAMETERS (cm)	FS9600	57	43	—	34	58×48×28	2.1	2	4.2
	FS9601	42	39	—	30	39×28×16.5	2.1	1	2.5
	FS9611	78	70	—	61-93	82×17×63	12.9	1	14.4
	FS9613	68	51-63	—	98-118	82×17×63	12.9	1	14.4

Walking Aids Series

• Foshan Wheelchair



95

FS910L FS920L FS930L FS9301L FS929L FS939L FS9206L FS9306L(14) FS9306L FS948L FS949L FS9203L



FS9275L-1 FS927L FS925LC FS933LC FS9332L FS933L FS9272LW FS937L



96



FS9300L FS938L(12) FS928L FS9281L FS9280L FS936L FS9361L FS583 FS938L



MAIN PARAMETERS (cm)	MODEL NO.	HEIGHT (cm)	N.W. (kg)			
			PCS CN	G.W. (kg)		
	FS910L	71-94	94×15×28	0.35	7.9	20
	FS920L	71-94	65×15×26	0.36	7.7	20
	FS930L	71-94	65×15×27	0.35	7.6	20
	FS9301L	71-94	65×15×28	0.35	7.6	20
	FS929L	71-94	65×15×26	0.36	7.8	20
	FS939L	67-90	63×15×28	0.34	7.3	20
	FS9206L	71-94	65×15×26	0.33	7	20
	FS9306L(14)	69-92	63×13×26	0.29	6.3	20
	FS9306L	69-92	63×13×26	0.29	6.3	20
	FS948L	76-86	61×16×26	0.37	8	20
	FS949L	76-86	61×16×26	0.37	8	20
	FS9203L	71-94	69×16×28	0.67	8	20
	FS9300L	75-98	74×33×15	0.4	8.8	20
	FS938L(12)	77-100	79×35×18	0.32	7.5	20
	FS928L	77-100	79×29×18	0.35	7.9	20
	FS9281L	75-98	69×35×18	0.33	7.5	20
	FS9280L	75-98	68×35×18	0.33	7.4	20
	FS936L	120	35×12×30	0.22	5	20
	FS9361L	120	38×12×30	0.3	6.4	20
	FS583	81	84×35×25	0.33	10.9	30
	FS938L	77-100	79×35×18	0.32	7.5	20
	FS9275L-1	83-93	72×16×26	0.37	8.3	20
	FS927L	80-90	61×16×26	0.37	8.3	20
	FS925LC	108-146	124×38×25	1.03	11.7	10
	FS933LC	102-121	70×33×32	0.58	12.8	20
	FS9332L	95-126	95×29×31	0.57	12.1	20
	FS933L	95-126	95×29×31	0.57	12.1	20
	FS9272LW	100-110	64×32×40	0.46	10.4	20
	FS937L	96-119	100×30×30	0.53	11.5	20
	FS937L	86-117	86×26×29	0.52	11.3	20
	FS925LB	112-143	109×32×27	0.72	15.6	20
	FS925LD	113-145	44×24×97	1	11.6	10

Walking Aids Series

• Foshan Wheelchair



97



98



FS9207LW-3

FS9207LW-2

FS9207LW-1

FS9304L

FS9381L

FS9382L



FS9242L

FS9415L

MAIN PARAMETERS (cm)	MODEL NO.	HIGHT (cm)	(cm)	N.W. (kg)	G.W. (kg)	PCS /CN
	FS9209L	72-95	67×16×28	0.37	8.2	20
	FS9391L	67-90	65×15×28	0.34	7.3	20
	FS9301L-1	71-94	61×35×29	0.39	8.6	20
	FS9301L-2	71-94	64×33×30	0.46	10.2	20
	FS9301L-3	71-94	66×15×28	0.38	8.4	20
	FS9301L-4	71-94	66×20×29	0.43	9.2	20
	FS9301L-5	71-94	64×30×43	0.46	10.2	20
	FS9208L	71-94	66×23×30	10.4	11.2	20
	FS9207LW-3	71-94	64×40×14	0.44	5.13	10
	FS9207LW-2	71-94	66×31×26	0.5	11.5	20
	FS9207LW-1	71-94	65×31×21	0.54	6.5	10
	FS9304L	94-117	117×24×42	0.79	17	20
	FS9381L	77-100	79×29×29	0.32	7.5	20
	FS9382L	77-100	79×40×28	0.51	11	20
	FS9411L	75-95	61×34×34	0.83	4.4	4
	FS9412L	70-90	116×29×43	0.81	9.6	10
	FS9413L	72-92	65×27×33	0.8	4.3	4
	FS9414L	71-94	78×41×35	0.63	7.5	10
	FS9243	66-89	71×22×45	0.91	8.4	8
	FS9221	73.5-96	78×29×41	0.83	9.5	10
	FS9242L	85-100	86×30×40	0.83	9.6	10
	FS9415L	75-95	61×34×34	0.8	4.3	4

Walking Aids Series

• Foshan Wheelchair



99



100



MAIN PARAMETERS (cm)	MODEL NO.	HIGHT (cm)		N.W. (kg)	G.W. (kg)	PCS CN
	FS926	74-97	78×29×36	0.76	8.7	10
	FS932	74-97	72×26×35	0.67	7.7	10
	FS942S	74-97	78×29×43	0.85	9.7	10
	FS945	73-88	78×29×43	0.78	8.8	10
	FS946S	74-97	78×29×43	0.85	9.5	10
	FS941	72-95	81×36×45	0.9	10.7	10
	FS9321	67-92	72×26×35	0.67	7.7	10
	FS9322	72-97	78×26×43	0.81	9.3	10
	FS924	65-88	72×30×34	0.8	9.2	10
	FS921L	74-97	83×41×46	0.74	9.4	10
	FS931L	74-97	83×41×46	0.74	9.4	10
	FS946A	74-97	78×29×43	0.85	9.5	10
	FS921	72-95	81×36×45	0.87	10.4	10
	FS931	72-95	81×36×45	0.9	10.6	10
	FS934	65-88	72×30×34	0.78	9.1	10
	FS947	72-95	78×29×36	0.77	8.8	10
	FS9213L	76-96	45×33×50	0.5	5.8	10
	FS924L	65-88	75×33×38	0.65	8	10
	FS9212L	70-88	78×29×36	0.5	6	10
	FS922	74-97	78×29×41	0.83	9.5	10
	FS911L	55	91×29×52	1.46	16.3	10

Walking Aids Series

Foshan Wheelchair



FS925L(L/M/S)



FS935(L/M/S)



FS943L

FS911L

FS945L

FS9113L

101

102



FS9204L

FS935S

FS9251L

FS9256L


FS9252L



FS940L

FS9112L

MAIN PARAMETERS (cm)	MODEL NO.	HIGHT (cm)	(cm)	N.W. (kg)	G.W. (kg)	PCS /CN
	FS925L(L)	135-155	140×33×28	0.88	19.5	20
	FS925L(M)	115-135	120×33×28	0.83	18.2	20
	FS925L(S)	95-115	101×33×28	0.73	16.5	20
	FS935(L)	128-148	135×34×28	0.85	18.9	20
	FS935(M)	113-133	120×34×28	0.75	16.5	20
	FS935(S)	98-118	106×34×28	0.7	15.3	20
	FS9204L	105-128	78×22×25	1.12	5.1	4
	FS935S(L)	119-139	130×20×30	1.26	16.5	12
	FS935S(M)	112-132	119×20×30	1.22	15.8	12
	FS935S(S)	95-115	101×20×30	1.2	15.3	12
	FS9251L	115-135	120×33×30	0.81	17.5	20
	FS9256(L)	119-139	125×33×29	1.5	31.8	20
	FS9256(M)	112-132	119×33×29	1.4	28.8	20
	FS9256(S)	95-115	85×29×40	0.87	9.7	10
	FS9252L	100-130	85×36×29	0.82	17.7	20
	FS911L	84	85×29×40	0.87	9.7	10
	FS943L	84	85×22×40	0.9	10	10
	FS945L	84	85×29×40	0.87	10	10
	FS9113L	84	85×29×40	0.87	9.7	10
	FS940L	80	81×20×41	0.85	9.3	10
	FS9112L	84	58×24×62	1.1	1.2	10




Hospital Bed Series

103



FS3241WGZF7

7 functions electric inclinable bed
With padel control system
Flat metal bed panel
Elevating abs guardrail with control
Powder coating steel frame
Removable abs head and foot board
Central control locking castor
Packing size: 219×109×63
Over all size: 218×108



FS3240WZF8

5 functions electric inclinable bed with
Padel control system
Flat metal bed panel with x-ray examination tray
Elevating abs guardrail
Detachable head and foot board
Coating steel frame
Central control locking castor
Weight measuring system, quick lie flat system
Packing size: 229×109×70
Over all size: 220×108



104



FS3238WGZF4-53W


5 functions electric inclinable bed with
Commode form
Net steel bed panel with mattress and
Commode bucket
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size: 200×102×56
Over all size: 207×96





FS3242WGZF8

5 functions electric inclinable bed
Flat metal bed panel
Elevating abs guardrail with control board
Powder coating steel frame
Abs head and foot board with control system
Central control locking castor
Weight measuring system
Packing size: 217×104×50
Over all size: 220×108



FS3133WBF6

3 functions electric bed with manual option
Detachable abs bed panel
Detachable full protection steel guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size: 205×96×43
Over all size: 212×102



FS3239WZF4

5 functions electric inclinable bed
Flat abs bed panel
Elevating abs guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Central control locking castor
Packing size: 224×104×49
Over all size: 228×108




FS3230WZ

3 functions electric bed
Net metal bed panel
Foldable aluminum guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Central control locking castor
Optional accessories: bed table
Packing size: 210×102×51
Over all size: 212×96



FS3238WGZF4

5 functions electric inclinable bed
Flat metal bed panel
Elevating guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Central control locking castor
Packing size: 217×103×45
Over all size: 210×108



FS3238W

5 functions electric inclinable bed
Net metal bed panel
Foldable aluminum guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size: 209×102×49
Over all size: 212×96



Hospital Bed Series

• Foshan Wheelchair



FS3233WG

3 functions electric bed with ultra-low position
Flat metal bed panel
Foldable aluminum guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size:219×102×32
Over all size:222×96

Max.150kg 22~64cm 0~80° 0~40°



FS3222WB

2 functions electric bed
Abs bed panel
Elevating abs guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size:209×95×47
Over all size:214×96

Max.175kg 45cm 0~80° 0~40°



FS3238WGZ4-54W

5 functions electric inclinable bed
Bed panel steel
Abs guardrail
Bed frame steel
Abs head and foot board
Pedal locking castor
Packing size:209×90×33
Over all size:210×100

Max.175kg 0~12° 52cm
0~70° 0~32° 0~58°
-40°~0°~40°



FS3236W

3 functions electric bed
Net metal bed panel
Foldable aluminum guardrail
Powder coating mild steel frame
Detachable abs head and foot board with Crashworthy wheels
Pedal locking castor
Packing size:218×102×49
Over all size:218×96

Max.175kg 45~78cm 0~80° 0~40°



FS3230W

3 functions electric bed
Net metal bed panel
Foldable aluminum guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size:209×102×49
Over all size:212×96


Max.175kg 45~78cm 0~80° 0~40°



FS3220W

2 functions electric bed
Net metal bed panel
Foldable aluminum guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size:209×102×28
Over all size:212×96

Max.175kg 45cm 0~80° 0~40°



Hospital Bed Series

109



Max.175kg 26~60cm 0~80° 0~40°



AU3230WMGZ

3 functions electric bed with linak Actuator
Flat metal bed panel
Foldable aluminum guardrail
Powder coating steel frame
Wooden head and foot board
Central control locking castor
Over all size:208×96
Linak actuator system



Max.175kg 32~72cm 0~80° 0~40°

FS3244WM

3 functions electric bed
Wooden bed panel
Elevating wooden guardrail
Powder coating mild steel frame
Wooden head and foot board
Pedal locking castor
Packing size:209×96×44
Over all size:209×103

110



Max.175kg 32~72cm 0~80° 0~40°

FS3246WM

3 functions electric bed
Wooden bed panel
Elevating wooden guardrail
Powder coating mild steel frame
Wooden head and foot board
Pedal locking castor
Packing size:209×96×30
Over all size:217×103



Max.175kg 38-70cm 0-76° 0-45°

111

FS3230WM

Electrical 3 func bed
Net metal panel
Foldable guard rails
Wooden head and foot board
Pedal locking castors



Max.175kg 46cm 0-75°

FS3210M-1

Electrical homcare bed
One function (backrest elevating)
Wooden frame
Wooden head and foot board
Wooden guard rails



Max.175kg 46cm 0-80° 0-40°

FS3220WM

Electrical 2 func bed
Net bed panel
Foldable guard rails
Wooden head and foot board
Pedal locking castors



Max.175kg 7-64cm 0-70° 0-20°

FS3237WMF4

Powder coating steel frame
Electric 3 funcs bed
No guard rails
Pedal locking castors
Low bed position

112



Max.175kg 45~78cm 0~80° 0~40°

113

FS3230WMGF3

3 functions electric bed
Flat metal bed panel
Full protection foldable aluminum
Guardrail
Wooden head and foot board
Powder coating steel frame
Pedal locking castor
Packing size:209×96×49
Over all size:208×96



Max.175kg 43~68cm 0~80° 0~40°

FS3231WM

3 functions electric bed with double
Motors lifting system
Net metal bed panel
Foldable aluminum guardrail
Powder coating steel frame
Wooden head and foot board
Pedal locking castor
Optional accessory: monkey pole
Packing size:215×101×34
Over all size:223×98


114



Max.175kg 45cm 0~80° 0~40°

FS3120WMB

2 functions electric bed with manual
Option
Abs bed panel
Foldable aluminum guardrail
Powder coating mild steel frame
Wooden head and foot board
Pedal locking castor
Packing size:204×102×33
Over all size:206×104



Max.175kg 45cm 0~80° 0~40°

FS3220MF2

2 functions electric bed,
Net metal bed panel
Detachable guardrail
Powder coating steel frame
Wooden head and foot board
Packing size:209×96×28
Over all size:208×97



FS3030W(FS3031W)

3 functions manual bed
Net metal panel
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size 209×92×25(209×102×41)
Over all size 212×96

Max.125kg 30~66cm (37~75cm) 0~80° 0~40°



FS3020W

2 fuctions manual bed
Net metal panel
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size 209×102×28
Over all size 212×96

Max.125kg 45cm 0~70° 0~30°



FS3026WG

2 fuctions manual bed
Flat metal panel
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size 224×102×32
Over all size 218×96

Max.125kg 48cm 0~75° 0~43°



FS3050W

5 fuctions manual bed with commode form
Turn over fuction on back section
Net steel bed panel with mattress and
Commode bucket
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size 211×98×55
Over all size 214×97

Max.125kg 51cm 0~75°
0~43° -40°~0°~40°



FS3010W

1 function manual bed
Net metal panel
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size 209×102×28
Over all size 212×96

Max.125kg 42cm 0~70°



FS3012W

1 fuction manual bed
Net metal panel
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size 209×102×28
Over all size 212×96

Max.125kg 41cm 0~80°



117

FS3131 Fysiocouch

Length: 200cm
Dynamic max load: 200kg
Static max load: 250kg
Pillow configurations
50cm-3-section; 70cm-3-section-
90cm-3-section; 70cm-4-section-
70cm-5-section; 70cm-2-section
(physiotherapy);
Adjustable height from 46cm to
95cm raised by actuator system
Five main pillows supported and
Angle-adjusted by gas springs
Armrests height and direction
Adjustable (optional parts)
Wheel lifting (central wheel lift,
Dual wheel lift)
Packing size: 174×75×60
Over all size: 170×71




FS3132 Fysiocouch

Length: 200cm
Dynamic max load: 200kg
Static max load: 250kg
Adjustable height from 46cm to
95cm raised by actuator system
Three main pillows supported and
Angle-adjusted by gas springs
Armrests height and direction
Adjustable (optional parts)
Wheel lifting (central wheel lift,
Dual wheel lift)
Packing size: 174×75×60
Overall size: 170×71



FS3020WG

2 functions manual bed
Flat metal
Foldable aluminum guardrail
Powder coating steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size 209×102×28
Over all size 212×96



FS3028WGF4

2 functions manual bed
Elevating abs guardrail
Powder coating mild steel frame
Detachable abs head and foot board
Pedal locking castor
Packing size:209×95×47
Over all size:214×96



FS3021SF2 Manual Stainless Steel Bed (FS3021F2 Manual Steel Bed)

2 functions manual bed
Detachable guardrail
Stainless steel frame
Stainless steel head and foot board
Over all size 205×93



118



Rehabilitation Aid Series



FS5771

- › Manual patient transfer chair, steel frame,
- › Solid castor,
- › Pedal locking rear wheel,
- › Fixed armrest,
- › Fixed footrest,
- › Commode seat,
- › Seat height manual adjusted

FS5772

- › Electric patient transfer chair,
- › Steel frame,
- › Solid castor,
- › Pedal locking rear wheel,
- › Fixed armrest,
- › Fixed footrest,
- › Commode seat,
- › Seat height electrical adjusted



AU5777

Electric patient hoist, powder coating steel frame
Maximum loading capacity: 150kg
Maximum height of lifting stand: 112cm
Width between legs: from 61cm to 92 cm



FS560

L: 88cm W:40cm H:73-125cm



FS5623

L: 90cm W: 45cm H: 72-113cm



FS5721

L: 77cm W: 38cm H: 72-112cm



FS572

L: 77cm W: 40cm H: 72-113cm

120



Rehabilitation Aid Series



FS3080W

Size: 187cm x 62cm x 65cm

FS5605S

Size: 180cm x 45cm x 174cm



FS5601S

Size: 83cm x 46cm x 64.5cm



FS5610S

Size: 190cm X 62.5cm X 99cm



FS5609S

Size: 72x45x81.5cm



FS5607S

Size: 72x48x81.5cm

Rehabilitation Aid Series

• Foshan Wheelchair



FS578B(1)

Size: 50cm×50cm×78cm



FS578B(2)

Size: 50cm×50cm×78cm



FS578B(3)

Size: 50cm×50cm×78cm



FS578B(4)

Size: 50cm×50cm×78cm



FS578

Size: 45cm X 47cm X 78cm



FS578H2



FS573

Size: L:90-118cm W: 33.5cm



FS575

Size: L:90 W: 30cm




FS590

Size: L:89cm W: 25cm



FS532

Mattress Size: 195 × 86 × 8cm



FS533

Mattress Size: 195 × 86 × 8cm



FS534

Mattress Size: 195 × 86 × 5cm



FS535

Mattress Size: 195 × 86 × 8cm



FS531



FS5311

123

124

Rehabilitation Aid Series

• Foshan Wheelchair



FS577

FS518

FS514

FS519

H:116cm-206cm

FS515

H:133cm-223cm

FS516

H:120cm-200cm

125



FS591-1

L×W: 31.5x18.5cm
weight capacity: 3kg

FS564

FS579



FS525G

FS525

FS525P



FS5606

FS565

FS512

FS566



FS507

FS561

FS505

FS506

126

Rehabilitation Aid Series

• Foshan Wheelchair



FS667B



FS8141



FS8142



D115(01)



D115(02)



D115(03)



FS661B



FS666B-100/50



FS664B



FS662B



FS663B




FS665B




FS7911S(12", 16", 24")




FS118900



FS7901B L300/L400/L500



FS7905SF1



PLASTIC ADAPTER




FS7904B

127

128

Spare Parts

• Foshan Wheelchair



Colour Sample

SEAT & BACK



Spare Parts



1XX: Spray-painting series 2XX: Powder coating series 3XX: Aluminum anodized series 4XX: Chromed series

(A6) 6"x1.5" Solid Castor
(A7) 5" Pvc Castor
(A8) 6" Pvc Castor
(B6) 8" Pvc Wh Eel Set
(B7) 8" Pneumatic Wheel Set
(B17) 5" Solid Wheelset
(C38) 24" Solid Spoke Wheel Set
(C41) 24" Solid Mag Wheel Set
(C42) 24" Pu Mag Wheel Set

(A26) 12" Pneumatic Mag Wheel Set
(A29) 16" Pneumatic Spoke Wheel Set
(A30) 12.5" 5-spoke Rubber Foam Rear Wheel
(A31) 12.5" 10-spoke Rubber Foam Rearwheel
(A32) 12" Rubber Foam Mag Wheel
(A33) 14" Pneumatic Rear Wheel
(C33) 16" Pu Wheelset
(C36) 24" Pneumatic Spokewheelset
(C37) 24" Pneumatic Mag Wheel Set
(C38) 24" Solid Spoke Wheel Set
(C41) 24" Solid Mag Wheel Set
(C42) 24" Pu Mag Wheel Set

(C43) 16" 3-spoke Pu Wheel Set
(C44) 24" 3-spoke Pu Wheel Set For Electric Wheelchair
(C45) Pu Rear Wheel With Drum Brake(20" , 22" , 24")
(C46) 24" Pneumatic Wheel Set For Electric Wheelchair
(E91) Steel Hub
(E96) Drum Brake Wheel Hub
(E97, 98)quick Release Hub
(E99) (I) 65 Drum Brake